

## Four Point Deep Flex Fixture

Four point flex fixture with three point option. Fixture is designed for testing high strength glass panels which exhibit tremendous deflection prior to break.

- Upper Anvils Articulate
- One Lower Anvil Articulates, One is Fixed.

## Specifications

Max Load	2 kN
Max Specimen Width	250mm/10in
Anvil Diameter	6 mm
Anvil Height	80mm
Mechanical Connections	Type Dm Lower, Type Om Upper

Note:  
Includes Alignment Plates for 20/40mm, 40/80mm, and 50/100mm Spans.

